

Abbott

Blue Laser 488nm

PMTs	Mirror	Filter	Detector	Fluorochromes
A	755 LP	780/60	B780	PE-Cy7
B	690 LP	710/50	B710	PerCP-Cy55
C	670 LP	675/20	B665	PE-Cy5
D	600 LP	610/20	B610	PE-Texas Red
E	550 LP	575/26	B575	PE
F	505 LP	525/50	B525	FITC
G	None	488/10	SSC	SSC
H				

Violet Laser 405nm

PMTs	Mirror	Filter	Detector	Fluorochromes
A	750 LP	780/60	V780	BV786
B	690 LP	710/50	V710	BV711
C	630 LP	660/60	V660	BV650
D	595 LP	610/20	V610	BV605
E	550 LP	585/42	V585	Qdot 585
F	535 LP	545/30	V545	Pacific Orange
G	505 LP	525/50	V525	BV510, Live dead Aqua
H	None	450/50	V450	Pacific Blue, BV421, DAPI

Red Laser 640nm

PMTs	Mirror	Filter	Detector	Fluorochromes
A	735LP	780/60	R780	APC-Cy7, APC-H7, APC-Fire750
B	685LP	710/40	R710	APC-Cy5.5, Alexa700, APC-R700
C	None	660/40	R660	APC, Alexa647, e660